

PLANETARY BALL MILL

LM PBM

For rapid fine crushing of soft, hard, brittle and fibrous material to end fineness $<1\mu\text{m}$

- Quick and easy to clean
- Rapid fine grinding
- Easy exchange of grinding jars and balls
- Grinding jars and balls made from a wide range of materials available
- Grinding jar volume from 50ml up to 500ml
- Can be expanded from 1 to 4 grinding jar positions
- CE-certified
- End fineness $<1\mu\text{m}$
- Programmable control



LAARMANN® Planetary Ball Mills for fine grinding of soft, hard to brittle or fibrous materials

The LAARMANN® Planetary Ball Mills are used for fine grinding of soft, hard to brittle or fibrous materials.

Dry and wet grindings are possible. They support the daily sample preparation for laboratory and development usage.

Working principle

LAARMANN® Planetary Ball Mills consist of several cylindrical grinding jars (positioned on the sun wheel as shown on the figure) which are filled with loose grinding balls. Two superimposed rotational movements move the grinding jars.

Like in a planetary system the grinding jar rotates on an orbit around the centre. This rotational movement is the self-rotation of the grinding container superimposed. The resulting centrifugal and acting acceleration forces lead to strong grinding effects.

Furthermore there are forces working according to the Coriolis acceleration. The result is an intensive grinding effect between the grinding balls and the sample.

There are different rotational ratios. With a rotation ratio of 1:-2 the grinding jar rotates twice during a sun wheel turn. The minus of this case indicate the opposite rotation direction.

Depending on the speed ratio different movement patterns of the grinding balls / media can be achieved. It can be achieved that the grinding media are crossing the grinding jar and loosen from the wall.

At hitting the wall of the grinding jar the sample will be stressed. At a different motion pattern the grinding balls roll over the sample and stress the ground material.

Operation

LAARMANN® Planetary Ball Mills enable the convenient programming of the following grinding parameters:


- Grinding duration in hours and minutes
- Speed
- Speed ratio (speed sun wheel : grinding jar)



PBM Planetary ball mill closed



Counter weight with 1 jar. can be expanded for 1-4 positions



Movement planetary mill

Features and benefits

- High efficient fine grinding up to end fineness $<1\mu\text{m}$
- Different speed ratios available
- Grinding jars from 50ml to 500 ml in different materials
- Suitable for long-term trials and continuous use
- Automatic direction reversal to avoid agglomerations
- Reproducible results due to programmable grinding parameters
- Can be expanded from 4 to 1 grinding jar positions
- CE - certified



Working principle

The selection of the right grinding jar and the correct filling level has a big impact on the grinding result. Depending on the application you have to select the correct material and amount/volume for the grinding jar and the grinding balls.

A jar filling should consist of about 1/3 sample and 1/3 ball charge. The remaining third is the free jar volume that is necessary for the movement of the balls. The following table provides recommendations.

We recommend to select always grinding jars and balls build from the same material

| Nominal volume | Sample amount | Recommended ball charge | | | | | | |
|----------------|---------------|-------------------------|-----------|----------|---------|---------|-------|------|
| | | Max feed size | 5mm | 10mm | 15mm | 20mm | 30mm | 40mm |
| 125 ml | 15-20 ml | $< 4 \text{ mm}$ | 500 pcs. | 30 pcs. | 18 pcs. | 7 pcs. | | |
| 250 ml | 25-120 ml | $< 6 \text{ mm}$ | 1200 pcs. | 50 pcs. | 45 pcs. | 15 pcs. | 6pcs. | |
| 500 ml | 75-220 ml | $< 10 \text{ mm}$ | 2000 pcs. | 100 pcs. | 70 pcs. | 25 pcs. | 8pcs. | 4pcs |



Grinding jars available with different liner:

- Stainless steel
- Hardened steel
- Tungsten carbide
- Agate
- Sintered corundum
- Zirconium oxide
- Pu-coated
- PTFE-coated

Features of grinding jars

The grinding jars are built from one block of steel or built with a stainless steel protective jacket with liner from the above materials.

Unique advantages of LAARMANN® grinding jars:

- Safe according to stainless steel protective jacket
- Easy opening according to gap between lid and jar
- Self-centring base of grinding jar
- Gas-tight and dust-proof sealed by o-ring
- Labelling field (e.g. for sample information)



500ml zirconium oxide




100ml agate



Grinding balls available in materials such as tungsten carbide, agate, zirconium oxide and stainless steel

Applications

- Floors
- Wood fibres
- Plant materials
- Seeds
- Tobacco
- Betonite
- Concrete
- Gypsum
- Sand
- Stone
- Cement clinker
- Hair
- Bones
- Tissue
- Carbon fibres
- Paints and lacquers
- Catalysts
- Plastics
- Pigments
- Polymers
- Cellulose
- Glass
- Hydroxylapatite
- Kaolin
- Ceramic oxides
- Quarz
- Clay minerals
- Ores
- Semi-precious stones
- Cole
- Coke
- Alloys
- Metal oxides
- Quarz
- Slags
- Electronic scrap
- Sludges
- Organic and unorganic waste



TECHNICAL DATA

| | |
|-------------------------|--|
| Electrical requirements | 200-240 Volt 50/60 Hz, others available on request |
| Motor power | 750 Watt |

TRANSPORT DATA

| | |
|------------------|-----------------------|
| Dimensions wxdxh | 759 x 479 x 548 mm |
| Weight | 95 kgs (only machine) |

PERFORMANCE

| | |
|-------------------------------------|---|
| Working principle | impact, friction |
| Feed size maximum | 10 mm |
| Number of grinding jars | 4/3/2/1 |
| Maximum volume of each grinding jar | 500 ml |
| Minimum volume of each grinding jar | 50 ml |
| Maximum end fineness | <1 μ m; <0,1 μ m for colloidal grinding |
| Adjustment of grinding duration | Digital adjustable |
| Start / Stop function | 1 |
| Setting of grinding time | Digital, 00:00 to 99:59 |
| Interval operation | Yes, with direction reversal |
| Pause time | 00:00 to 99:59 |
| Speed ratio | 1:-2 |
| Sun wheel speed | 50-650 min ⁻¹ |
| Effective sun wheel diameter | 360 mm |